

COOPERS®
Guardian® 6 in 1
Vaccine for Sheep and Lambs



For immunity against

- **Cheesy Gland** (*Corynebacterium pseudotuberculosis*)
- **Pulpy Kidney** (*Clostridium perfringens* Type D)
- **Blackleg** (*Clostridium chauvoei*)
- **Malignant Oedema** (*Clostridium septicum*)
- **Black Disease** (*Clostridium novyi* Type B)
- **Tetanus** (*Clostridium tetani*)

in sheep and lambs.

250 DOSES

PM 49
C/C 90244

FOR ANIMAL TREATMENT ONLY

COOPERS®
Guardian® 6 in 1
Vaccine for Sheep and Lambs

ACTIVE CONSTITUENTS: *Corynebacterium pseudotuberculosis* toxoid and inactivated cells.
Clostridium perfringens Type D toxoid. *Clostridium chauvoei* toxoid and inactivated cells.
Clostridium septicum toxoid. *Clostridium novyi* Type B toxoid and inactivated cells.
Clostridium tetani toxoid.

Contains 0.015 mg/mL thiomersal as preservative.

For immunity against Cheesy Gland, Pulpy Kidney, Blackleg, Malignant Oedema, Black Disease and Tetanus in sheep and lambs.



COOPERS



500mL

COOPERS®
Guardian® 6 in 1
Vaccine for Sheep and Lambs



SPECIALIST ADVICE IN EMERGENCY ONLY
Schering-Plough Pty Ltd
1 800 226 511
ALL HOURS - AUSTRALIA WIDE

® Schering-Plough Pty Ltd Registered Trademark
© Schering-Plough Pty Ltd 2001

250 DOSES

AMCOR CARTONS
AUSTRALIA
12
PD 101063
ABCDEFGHIJKLMN
OPQRSTUWXYZ

COOPERS®
Guardian® 6 in 1
Vaccine for Sheep and Lambs



GLOSS FREE AREA

READ THE ENCLOSED LEAFLET BEFORE USING THIS PRODUCT.


DIRECTIONS FOR USE

Shake contents before use. USE ALL PRODUCT WITHIN 24 HOURS OF OPENING.

Do not expose directly to heat and light during use. The use of an insulated container is recommended during the day to protect the vaccine. If it is to be hung beside the race, ensure it is on the shady side.

Dose: Sheep and Lambs: 2mL by subcutaneous injection.

1 Sensitising Dose **2 2nd Dose (1st Booster): After 4-6 weeks** **3 Annual Booster**



Programme: Sensitising Dose, followed 4-6 weeks later by a Booster Dose. Then a Booster Dose annually.

THE VACCINATION PROCEDURE

CAUTION: AVOID CARCASS DAMAGE

1. Sterilise all injection apparatus by boiling for at least ten minutes before use. Avoid use of strong disinfectants on apparatus.
2. Maintain maximum cleanliness at all times.

3. Keep needles sharp and clean. Change needles frequently (every 20-30 animals). Use short needles. Recommended maximum sizes 15mm x 18 gauge for sheep and lambs.

4. Check that the equipment is functioning properly and that it is delivering the correct dose.
5. As far as possible avoid injection of animals during wet weather or under dusty conditions.
6. Inject only under the skin and not into the muscle, by raising a fold of skin and injecting the dose under the raised skin - **NOT INTO IT.**
7. Sheep and lambs should be injected high on the side of the neck, behind the ear.

WITHHOLDING PERIOD NIL

DISPOSE of empty container by wrapping in paper and putting in garbage.

STORE at 2 to 8°C (REFRIGERATE, DO NOT FREEZE)

Keep out of reach of children.

- Protect from Light
- APVMA Approval No 46578/1199

EXP. REFER TO BATCH AND EXPIRY ON TOP PANEL



COOPERS ANIMAL HEALTH
A division of Schering-Plough Pty Ltd
11 Gibbon Road, Baulkham Hills NSW 2153

500mL

